

# GL CATALOG

## 2024



# UPPER EXTREMITY DEVICES

## GL Finger Device

Each complete JAS GL Finger Device requires: 1 – Finger Frame, 1 – Finger Hand Cuff  
 Each complete JAS GL Thumb Device requires: 1 – Thumb Frame, 1 – Thumb Hand Cuff



Description	Product Code	Finger Circumference (D) in. (cm)
GL Finger Frame, Medium	GL50M-000-00	Up to 2 1/4 ( 5.7)
GL Finger Frame, Large	GL50N-000-00	2 1/4 (5.7) and up
GL Thumb Frame, Medium	GL500-00M-T0	Up to 2 1/4 ( 5.7)
GL Thumb Frame, Large	GL500-00L-T0	2 1/4 (5.7) and up

ROM: 100° Flexion – 10° Hyperextension

**GL Finger Hand Cuffs** (See Upper Extremity Measurement Chart for measurement locations)

Hand Cuff Size	Product Code		Width of Hand (H) in. (cm)
	Left	Right	
Small	GL500-000-SL1	GL500-000-SR1	up to 3 (7.6)
Medium	GL500-000-ML1	GL500-000-MR1	3 – 3 1/2 (7.6 – 8.9)
Large	GL500-000-LL1	GL500-000-LR1	3 1/2 - 4 (8.9 – 10.2)
XLarge	GL500-000-XLL1	GL500-000-XLR1	4 (10.2) and up

**GL Thumb Hand Cuffs:** (No measurements required)

Description	Product Code	
	Left	Right
GL Thumb Hand Cuff	GL500-000-TL	GL500-000-TR

## GL Wrist Device

Each complete JAS GL Wrist Device requires: 1 – Frame, 1 – Forearm Cuff

Description	Product Code	
	Left	Right
GL Wrist Frame	GL400-000-L1	GL400-000-R1



**GL Wrist Forearm Cuffs<sup>1</sup>** (See Upper Extremity Measurement Chart for measurement locations) ROM: 10° Flexion – 90° Extension

Forearm Cuff Size	Product Code		Forearm Circumference (B) in. (cm)	Forearm Length (F) in. (cm)
	Left	Right		
Small, Short	GL740-SXN-FL	GL740-SXN-FR	7 -8 1/4 (17.8 – 21.0)	
Medium, Short	GL740-SN6-FL	GL740-SN6-FR	8 1/4 -10 1/8 (21.0 - 25.7)	
Large, Short	GL740-SW6-FL	GL740-SW6-FR	10 1/8 – 12 1/4 (25.7 – 31.1)	7 - 9 (17.8 – 22.9)
Xlarge, Short	GL740-SXW-FL	GL740-SXW-FR	12 1/4 - 14 1/8 (31.1 – 35.9)	
2Xlarge, Short	GL740-SXXW-FL	GL740-SXXW-FR	14 1/8 (35.9) and up	
Small, Long	GL740-MXN-FL	GL740-MXN-FR	7 -8 1/4 (17.8 – 21.0)	
Medium, Long	GL740-MN6-FL	GL740-MN6-FR	8 1/4 -10 1/8 (21.0 - 25.7)	
Large, Long	GL740-MW6-FL	GL740-MW6-FR	10 1/8 – 12 1/4 (25.7 – 31.1)	9 (22.9) and up
Xlarge, Long	GL740-MXW-FL	GL740-MXW-FR	12 1/4 - 14 1/8 (31.1 – 35.9)	
2Xlarge, Long	GL740-MXXW-FL	GL740-MXXW-FR	14 1/8 (35.9) and up	

1.) Includes: Hand Pad, Hand Strap

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# UPPER EXTREMITY DEVICES

## GL Pronation/Supination Device

Each complete JAS GL Pro/Sup Device requires: 1 – Frame, 1 – Bicep Cuff, 1 – Hand Cuff



ROM: 110° Pronation – 130° Supination

Description	Product Code		Forearm Length (G) in. (cm)
	Left	Right	
GL Pro/Sup Frame, Short	GL060-00M-L1	GL060-00M-R1	6 - 10 (15.2 – 25.4)
GL Pro/Sup Frame, Long	GL060-00L-L1	GL060-00L-R1	10 - 14 (25.4 -35.6)

### Pro/Sup Cuffs:

**GL Pro/Sup Bicep Cuffs<sup>2</sup>** (See Upper Extremity Measurement Chart for measurement locations)

Bicep Cuff Size	Product Code		Bicep Circumference (A) in. (cm)	Bicep Length (E) in. (cm)
	Left	Right		
XSmall	GL706-XSE-BL	GL706-XSE-BR	7 1/4 - 8 3/4 (18.4 – 22.2)	6 7/8 (17.5 ) and up
Small	GL706-S0E-BL	GL706-S0E-BR	8 3/4 – 10 (22.2 – 25.4)	
Medium	GL706-M0E-BL	GL706-M0E-BR	10 – 11 1/4 (25.4 – 28.6)	
Large	GL706-L0E-BL	GL706-L0E-BR	11 1/4 - 12 1/4 (28.6 – 31.1)	
XLarge	GL706-XLE-BL	GL706-XLE-BR	12 1/4 - 13 3/8 (31.1 - 34.0)	
2XLarge	GL706-XXLE-BL	GL706-XXLE-BR	13 3/8 - 15 1/2 (34.0 – 39.4)	
3XLarge	GL706-3XLE-BL	GL706-3XLE-BR	15 1/2 - 17 (39.4 – 43.2)	
4XLarge	GL706-4XLE-BL	GL706-4XLE-BR	17 – 17 3/4 (43.2 – 45.1)	
5XLarge	GL706-5XLE-BL	GL706-5XLE-BR	17 3/4 - 19 1/4 (45.1 – 48.9)	
6XLarge	GL706-6XLE-BL	GL706-6XLE-BR	19 1/4 – 21 (48.9 – 53.3)	

2.) Includes: Elbow Pad

**GL Pro/Sup Hand Cuffs** (See Upper Extremity Measurement Chart for measurement locations)

Hand Cuff Size	Product Code		Wrist Circumference (C) in. (cm)
	Left	Right	
1, Short	GL706-1PH-L1	GL706-1PH-R1	5 - 5 3/4 (12.7 - 14.6)
5, Short	GL706-5PH-L1	GL706-5PH-R1	5 3/4 - 7 (14.6 - 17.8)
8, Short	GL706-8PH-L1	GL706-8PH-R1	7 - 7 1/2 (17.8 - 19.1)
10, Short	GL706-10PH-L1	GL706-10PH-R1	7 1/2 - 8 1/4 (19.1 - 21.0)

**GL Pro/Sup Extended Hand Cuffs** (See Upper Extremity Measurement Chart for measurement locations)

Hand Cuff Size	Product Code		Wrist Circumference (C) in. (cm)
	Left	Right	
1, Long	GL706-1PH-LE	GL706-1PH-RE	5 - 5 3/4 (12.7 - 14.6)
5, Long	GL706-5PH-LE	GL706-5PH-RE	5 3/4 - 7 (14.6 - 17.8)
8, Long	GL706-8PH-LE	GL706-8PH-RE	7 - 7 1/2 (17.8 - 19.1)
10, Long	GL706-10PH-LE	GL706-10PH-RE	7 1/2 - 8 1/4 (19.1 - 21.0)



ROM: 110° Pronation – 130° Supination

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# UPPER EXTREMITY DEVICES

## GL Elbow Device

Each complete JAS GL Elbow Device requires: 1 – Frame, 1 – Bicep Cuff, 1 – Forearm Cuff



ROM: 167° Flexion – 15° Hyperextension

Description	Product Code	
	Left	Right
GL Elbow Frame	GL100-000-L1	GL100-000-R1

### GL Elbow Cuffs:

**GL Elbow Bicep Cuffs** (See Upper Extremity Measurement Chart for measurement locations)

Bicep Cuff Size	Product Code		Bicep Circumference (A) in. (cm)	Bicep Length (E) in. (cm)
	Left	Right		
XSmall	GL710-XSE-BL	GL710-XSE-BR	7 1/4 - 8 3/4 (18.4 – 22.2)	7 1/2 (19.0) and up
Small	GL710-S0E-BL	GL710-S0E-BR	8 3/4 – 10 (22.2 – 25.4)	
Medium	GL710-M0E-BL	GL710-M0E-BR	10 – 11 1/4 (25.4 – 28.6)	
Large	GL710-L0E-BL	GL710-L0E-BR	11 1/4 - 12 1/4 (28.6 – 31.1)	
XLarge	GL710-XLE-BL	GL710-XLE-BR	12 1/4 - 13 3/8 (31.1 - 34.0)	
2XLarge	GL710-XXLE-BL	GL710-XXLE-BR	13 3/8 - 15 1/2 (34.0 – 39.4)	
3XLarge	GL710-3XLE-BL	GL710-3XLE-BR	15 1/2 - 17 (39.4 – 43.2)	
4XLarge	GL710-4XLE-BL	GL710-4XLE-BR	17 – 17 3/4 (43.2 – 45.1)	
5XLarge	GL710-5XLE-BL	GL710-5XLE-BR	17 3/4 - 19 1/4 (45.1 – 48.9)	
6XLarge	GL710-6XLE-BL	GL710-6XLE-BR	19 1/4 – 21 (48.9 – 53.3)	

**GL Elbow Forearm Cuffs** (See Upper Extremity Measurement Chart for measurement locations)

Forearm Cuff Size	Product Code		Forearm Circumference (B) in. (cm)	Forearm Length (G) in. (cm)
	Left	Right		
XSmall, Short	GL710-XSS-FL	GL710-XSS-FR	7 - 8 (17.8 – 20.3)	7 - 9 1/2 (17.8 - 24.1)
Small, Short	GL710-S0S-FL	GL710-S0S-FR	8 - 10 (20.3 - 25.4)	
Medium, Short	GL710-M0S-FL	GL710-M0S-FR	10 - 12 (25.4 - 30.5)	
Large, Short	GL710-L0S-FL	GL710-L0S-FR	12 - 14 (30.5 – 35.6)	
XLarge, Short	GL710-XLS-FL	GL710-XLS-FR	14 (35.6) and up	
XSmall, Long	GL710-XSL-FL	GL710-XSL-FR	7 - 8 (17.8 – 20.3)	9 1/2 (24.1) and up
Small, Long	GL710-S0L-FL	GL710-S0L-FR	8 - 10 (20.3 - 25.4)	
Medium, Long	GL710-M0L-FL	GL710-M0L-FR	10 - 12 (25.4 - 30.5)	
Large, Long	GL710-L0L-FL	GL710-L0L-FR	12 - 14 (30.5 – 35.6)	
XLarge, Long	GL710-XLL-FL	GL710-XLL-FR	14 (35.6) and up	

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# UPPER EXTREMITY DEVICES

## GL Shoulder Device: Body Cuff

Each complete JAS Shoulder Body Cuff Device requires: 1 – Frame,  
1 – Body Cuff Cuff Set, 1 - Chest Cuff, 1 – Pelvic Cuff



ROM: 100° External Rotation –  
90° Internal Rotation, 90° Abduction

Description	Product Code		Length (E) in. (cm)
	Left	Right	
GL Shoulder Frame, Short, Body Cuff	GL070-00S-LB	GL070-00S-RB	7 1/2 – 8 3/4 (19 – 22.2)
GL Shoulder Frame, Long, Body Cuff	GL070-00L-LB	GL070-00L-RB	8 5/8 – 10 3/4 (21.9 – 27.3)

### GL Shoulder Cuffs:

#### GL Shoulder Body Cuff: Cuff Set<sup>3</sup>

Product Code	
Left	Right
GL707-B00-L0	GL707-B00-R0

3.) Includes: Bicep Cuff, Forearm Cuff, Hand Cuff

#### GL Shoulder Chest Cuffs (See Upper Extremity Measurement Chart for measurement locations)

Chest Cuff Size	Product Code		Chest Circumference (I) in. (cm)
	Left	Right	
Small	GL707-00S-CL	GL707-00S-CR	28 -34 3/4 (71.1 – 88.3)
Medium	GL707-00M-CL	GL707-00M-CR	34 3/4 – 41 1/2 (88.3 – 105.4)
Large	GL707-00L-CL	GL707-00L-CR	41 1/2 – 48 1/4 (105.4 – 122.6)
XLarge	GL707-0XL-CL	GL707-0XL-CR	48 1/4 (122.6) and up

#### GL Shoulder Pelvic Cuffs (See Upper Extremity Measurement Chart for measurement locations)

Pelvic Cuff Size	Product Code		Waist Circumference (J) in. (cm)
	Left	Right	
Small	GL707-00S-PL	GL707-00S-PR	24 -33 3/4 (61.0 – 85.7)
Medium	GL707-00M-PL	GL707-00M-PR	33 3/4 - 42 1/2 (85.7 – 108.0)
Large	GL707-00L-PL	GL707-00L-PR	42 1/2 - 51 1/4 (108.0 – 130.2)
XLarge	GL707-0XL-PL	GL707-0XL-PR	51 1/4 (130.2) and up

**Floor Stand Option:** Order a Floor Stand instead of a Pelvic Cuff.


Description	Product Code
Floor Stand	S07GL-003-01

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# UPPER EXTREMITY DEVICES

## Upper Extremity Measurement Chart – GL Devices

- A Circumference: Largest point of Bicep
- B Circumference: Largest point of Forearm
- C Circumference: 1 1/2" (3.8 cm) proximal to Ulnar Styloid
- D Circumference: Affected Finger proximal Phalanx
- E Length: Axilla to Medial Epicondyle
- F Length: Crease of Elbow to Ulnar Styloid
- G Length: Lateral Epicondyle to Ulnar Styloid
- H Length: Width of Hand across MP Joints
- I Circumference: Chest at Axilla
- J Circumference: Waist at Belly Button



**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# LOWER EXTREMITY DEVICES

## GL Knee Extension Device

Each complete JAS GL Knee Extension Device requires: 1 – Frame,  
1 – Thigh Cuff, 1 – Calf Cuff



ROM: 48° Flexion – 36° Hyperextension

Description	Product Code	
	Left	Right
GL Knee Extension Frame	GL200-000-L0	GL200-000-R0

### GL Knee Extension Cuffs:

**GL Knee Extension Thigh Cuffs** (See Lower Extremity Measurement Chart for measurement locations)

Thigh Cuff Size	Product Code		Thigh Circumference (A) in. (cm)	Thigh Length (C) in. (cm)
	Left	Right		
Small, Short	GL720-ESS-TL	GL720-ESS-TR	14 -17 (35.6 - 43.2)	9 - 12 (22.9 – 30.5)
Small, Long	GL720-ESL-TL	GL720-ESL-TR		
Medium, Short	GL720-EMS-TL	GL720-EMS-TR	17 - 20 (43.2 - 50.8)	9 - 12 (22.9 – 30.5)
Medium, Long	GL720-EML-TL	GL720-EML-TR		
Large, Short	GL720-ELS-TL	GL720-ELS-TR	20 - 23 (50.8 - 58.4)	9 - 12 (22.9 – 30.5)
Large, Long	GL720-ELL-TL	GL720-ELL-TR		
XLarge, Short	GL720-EXLS-TL	GL720-EXLS-TR	23 - 26 (58.4 - 66.0)	9 - 12 (22.9 – 30.5)
XLarge, Long	GL720-EXLL-TL	GL720-EXLL-TR		
2XLarge, Short	GL720-E2XLS-TL	GL720-E2XLS-TR	26 - 29 (66.0 - 73.7)	9 - 12 (22.9 – 30.5)
2XLarge, Long	GL720-E2XLL-TL	GL720-E2XLL-TR		
3XLarge, Short	GL720-E3XLS-TL	GL720-E3XLS-TR	29 - 33 (73.7 – 83.8)	9 - 12 (22.9 – 30.5)
3XLarge, Long	GL720-E3XLL-TL	GL720-E3XLL-TR		

**GL Knee Extension Calf Cuffs** (See Lower Extremity Measurement Chart for measurement locations)


Calf Cuff Size	Product Code		Calf Circumference (B) in. (cm)	Tibia Length (D) in. (cm)
	Left	Right		
Small, Short	GL720-ESS-CL	GL720-ESS-CR	11 -14 (27.9 - 35.6)	13 – 16 1/2 (30.5 – 41.9)
Small, Long	GL720-ESL-CL	GL720-ESL-CR		
Medium, Short	GL720-EMS-CL	GL720-EMS-CR	14 - 16 (35.6 – 40.6)	13 – 16 1/2 (30.5 – 41.9)
Medium, Long	GL720-EML-CL	GL720-EML-CR		
Large, Short	GL720-ELS-CL	GL720-ELS-CR	16 - 18 (40.6 – 45.7)	13 – 16 1/2 (30.5 – 41.9)
Large, Long	GL720-ELL-CL	GL720-ELL-CR		
XLarge, Short	GL720-EXLS-CL	GL720-EXLS-CR	18 - 20 (45.7 – 50.8)	13 – 16 1/2 (30.5 – 41.9)
XLarge, Long	GL720-EXLL-CL	GL720-EXLL-CR		

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# LOWER EXTREMITY DEVICES

## GL Knee Flexion Device

Each complete JAS GL Knee Flexion Device requires:  
1 – Frame, 1 – Thigh Cuff, 1 – Tibia Cuff



ROM: 50° Flexion – 148° Hyperextension

Description	Product Code
GL Knee Flexion Frame	GL20FX-000-01

### GL Knee Flexion Cuffs:

**GL Knee Flexion Thigh Cuffs** (See Lower Extremity Measurement Chart on for measurement locations)

Thigh Cuff Size	Product Code		Thigh Circumference (A) in. (cm)	Thigh Length (C) in. (cm)
	Left	Right		
Small, Short	GL720-FSS-TL	GL720-FSS-TR	14 -17 (35.6 - 43.2)	9 - 12 (22.9 – 30.5)
Small, Long	GL720-FSL-TL	GL720-FSL-TR		12 (30.5) and up
Medium, Short	GL720-FMS-TL	GL720-FMS-TR	17 - 20 (43.2 - 50.8)	9 - 12 (22.9 – 30.5)
Medium, Long	GL720-FML-TL	GL720-FML-TR		12 (30.5) and up
Large, Short	GL720-FLS-TL	GL720-FLS-TR	20 - 23 (50.8 - 58.4)	9 - 12 (22.9 – 30.5)
Large, Long	GL720-FLL-TL	GL720-FLL-TR		12 (30.5) and up
XLarge, Short	GL720-FXLS-TL	GL720-FXLS-TR	23 - 26 (58.4 - 66.0)	9 - 12 (22.9 – 30.5)
XLarge, Long	GL720-FXLL-TL	GL720-FXLL-TR		12 (30.5) and up
2XLarge, Short	GL720-F2XLS-TL	GL720-F2XLS-TR	26 - 29 (66.0 - 73.7)	9 - 12 (22.9 – 30.5)
2XLarge, Long	GL720-F2XLL-TL	GL720-F2XLL-TR		12 (30.5) and up
3XLarge, Short	GL720-F3XLS-TL	GL720-F3XLS-TR	29 - 33 (73.7 – 83.8)	9 - 12 (22.9 – 30.5)
3XLarge, Long	GL720-F3XLL-TL	GL720-F3XLL-TR		12 (30.5) and up

**GL Knee Flexion Tibia Cuffs** (See Lower Extremity Measurement Chart on for measurement locations)

Tibia Cuff Size	Product Code		Calf Circumference (B) in. (cm)	Tibia Length (D) in. (cm)
	Left	Right		
Small, Short	GL720-FSS-CL	GL720-FSS-CR	11 -14 (27.9 – 35.6)	13 – 16 1/2 (30.5 – 41.9)
Small, Long	GL720-FSL-CL	GL720-FSL-CR		16 1/2 (41.9) and up
Medium, Short	GL720-FMS-CL	GL720-FMS-CR	14 - 16 (35.6 – 40.6)	13 – 16 1/2 (30.5 – 41.9)
Medium, Long	GL720-FML-CL	GL720-FML-CR		16 1/2 (41.9) and up
Large, Short	GL720-FLS-CL	GL720-FLS-CR	16 - 18 (40.6 – 45.7)	13 – 16 1/2 (30.5 – 41.9)
Large, Long	GL720-FLL-CL	GL720-FLL-CR		16 1/2 (41.9) and up
XLarge, Short	GL720-FXLS-CL	GL720-FXLS-CR	18 - 20 (45.7 – 50.8)	13 – 16 1/2 (30.5 – 41.9)
XLarge, Long	GL720-FXLL-CL	GL720-FXLL-CR		16 1/2 (41.9) and up
2XLarge, Short	GL720-F2XLS-CL	GL720-F2XLS-CR	20 - 22 (50.8 – 55.9)	13 – 16 1/2 (30.5 – 41.9)
2XLarge, Long	GL720-F2XLL-CL	GL720-F2XLL-CR		16 1/2 (41.9) and up

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**



# LOWER EXTREMITY DEVICES

## GL Ankle Device

Each complete JAS GL Ankle Device requires: 1 – Frame, 1 – Calf Cuff



Description	Product Code	
	Left	Right
GL Ankle Frame	GL050-000-L1	GL050-000-R1

ROM: 40° Dorsiflexion – 52° Plantarflexion

**GL Ankle Cuffs<sup>5</sup>** (See Lower Extremity Measurement Chart for measurement locations.)

Calf Cuff Size	Product Code		Calf Circumference (B) in. (cm)	Tibia Length (D) in. (cm)
	Left	Right		
Small, Short	GL705-0SS-L1	GL705-0SS-R1	11 -13 (27.9 – 33.0)	11 -14 1/2 (27.9 – 33.8)
Small, Long	GL705-0SL-L1	GL705-0SL-R1		
Medium, Short	GL705-0MS-L1	GL705-0MS-R1	13 - 15 (33.0 – 38.1)	11 -14 1/2 (27.9 – 33.8)
Medium, Long	GL705-0ML-L1	GL705-0ML-R1		
Large, Short	GL705-0LS-L1	GL705-0LS-R1	15 - 17 (38.1 – 43.2)	11 -14 1/2 (27.9 – 33.8)
Large, Long	GL705-0LL-L1	GL705-0LL-R1		
XLarge, Short	GL705-XLS-L1	GL705-XLS-R1	17 - 20 (43.2 – 50.8)	11 -14 1/2 (27.9 – 33.8)
XLarge, Long	GL705-XLL-L1	GL705-XLL-R1		

5.) Includes: Heel Cuff, Foot Plate Pad, Foot Plate Strap

## GL Toe Device

Each complete JAS GL Toe Device requires: 1 – Frame, 1 – Foot Plate



ROM: 90° Extension – 45° Flexion

Description	Product Code
GL Toe Frame	GL030-000-01

**GL Toe Foot Plates<sup>5</sup>** (See Lower Extremity Measurement Chart for measurement locations.)

Foot Plate Size	Product Code		Foot Length (E) in. (cm)
	Left	Right	
Small	GL030-S00-L0	GL030-S00-R0	up to 9 1/2 (24.1)
Medium	GL030-M00-L0	GL030-M00-R0	9 1/2 – 10 1/2 (24.1 – 26.7)
Large	GL030-L00-L0	GL030-L00-R0	10 1/2 – 11 5/16 (26.7 – 28.7)
XLarge	GL030-XL0-L0	GL030-XL0-R0	11 5/16 – 12 5/16 (28.7 – 31.3)


5.) Includes: Upright

**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**

# LOWER EXTREMITY DEVICES

## Lower Extremity Measurement Chart – GL Devices

- A Circumference: 1 1/2" (3.8 cm) below Groin
- B Circumference: Largest point of Calf
- C Length: Groin to Knee Center (*not center of Knee Cap*)
- D Length: Knee Center (*not center of Knee Cap*) to Medial Malleolus
- E Length: Foot Length



**NOTE: MODIFICATIONS AND CUSTOMIZED SIZES OUTSIDE THE CATALOG SIZES MAY INCLUDE AN UPCHARGE**