

# **PEDIATRIC CATALOG**

## **2017**



# UPPER EXTREMITY DEVICES

## Pediatric Wrist Device

Each complete JAS Pediatric Wrist Device requires: 1 – Frame, 1 – Forearm Cuff



Description	Product Code
Pediatric Wrist Frame	0040P-000-00

### Pediatric Wrist Forearm Cuffs<sup>1</sup> (See Upper Extremity Measurement Chart on page 5 for measurement locations)

Forearm Cuff Size	Product Code	Forearm Circumference (B) in. (cm)	Forearm Length (F) cm (in.)
Pediatric	00740-P00-F0	3 1/4 – 4 1/4 (8.3 – 10.8)	Greater than 5 (13.3)
Small	00740-P0S-F0	4 1/4 – 5 3/4 (10.8 – 14.6)	Greater than 6 1/2 (15.2)
Medium	00740-P0M-F0	5 3/4 – 7 3/4 (14.6 – 20.0)	Greater than 6 1/2 (15.2)

1.) Includes: Hand Pad, Hand Strap

# UPPER EXTREMITY DEVICES

## Pediatric Pronation/Supination Device

Each complete JAS Pediatric Pro/Sup Device requires: 1 – Frame, 1 – Bicep Cuff, 1 – Hand Cuff (Included with Bicep Cuff)



Description	Product Code		Device Size
	Left	Right	
Pediatric Pro/Sup Frame	0006P-000-L0	0006P-000-R0	Forearm Length (G) is less than 8" (20.3 cm) and greater than or equal to 5 1/2" (14.0 cm).

### Pro/Sup Cuffs:

**Pediatric Pro/Sup Bicep Cuffs<sup>2</sup>** (See Upper Extremity Measurement Chart on page 5 for measurement locations)

Bicep Cuff Size	Product Code		Bicep Circumference (A) in. (cm)	Bicep Length (E) in. (cm)
	Left	Right		
Pediatric	00706-P00-BL	00706-P00-BR	5 3/4 – 7 1/4 (14.6 – 18.4)	
Xsmall	00706-PXS-BL	00706-PXS-BR	7 1/4 - 8 3/4 (18.4 – 22.2)	Greater than 5 1/2 (14.0)
Small	00706-P0S-BL	00706-P0S-BR	8 3/4 – 10 (22.2 – 25.4)	

2.) Includes: Hand Cuff

# UPPER EXTREMITY DEVICES

## Pediatric Elbow Device

Each complete JAS Pediatric Elbow Device requires: 1 – Frame, 1 – Bicep Cuff, 1 – Forearm Cuff



Description	Product Code
Pediatric Elbow Frame	0010P-000-00

### Pediatric Elbow Cuffs:

**Pediatric Elbow Bicep Cuffs** (See Upper Extremity Measurement Chart on page 5 for measurement locations)

Bicep Cuff Size	Product Code		Bicep Circumference (A) in. (cm)	Bicep Length (E) in. (cm)
	Left	Right		
Pediatric	00710-P00-BL	00710-P00-BR	4 3/4 – 7 1/4 (12.1 – 18.4)	6 1/2 (15.2) and up
Xsmall	00710-PXS-BL	00710-PXS-BR	7 1/4 - 8 3/4 (18.4 – 22.2)	7 (17.8) and up
Small	00710-P0S-BL	00710-P0S-BR	8 3/4 – 10 (22.2 – 25.4)	7 (17.8) and up


**Pediatric Elbow Forearm Cuffs** (See Upper Extremity Measurement Chart on page 5 for measurement locations)

Forearm Cuff Size	Product Code		Forearm Circumference (B) in. (cm)	Forearm Length (G) in. (cm)
	Left	Right		
Pediatric	00710-P00-FL	00710-P00-FR	4 1/2 - 5 3/4 (10.8 – 14.6)	6 1/2 (15.2) and up
Small	00710-P0S-FL	00710-P0S-FR	5 3/4 – 7 3/4 (14.6 – 19.7)	8 1/4 (21.0) and up
Medium	00710-P0M-FL	00710-P0M-FR	7 3/4 – 10 3/4 (19.7 – 27.3)	9 (22.9) and up

# UPPER EXTREMITY DEVICES

## Upper Extremity Measurement Chart – Pediatric Devices

- A Circumference: Largest point of Bicep
- B Circumference: Largest point of Forearm
- C Circumference: 1 1/2" (3.8 cm) proximal to Ulnar Styloid
- D Circumference: Affected Finger proximal Phalanx
- E Length: Axilla to Medial Epicondyle
- F Length: Crease of Elbow to Ulnar Styloid
- G Length: Lateral Epicondyle to Ulnar Styloid
- H Length: Width of Hand across MP Joints
- I Circumference: Chest at Axilla
- J Circumference: Waist at Belly Button



# LOWER EXTREMITY DEVICES

## Pediatric Knee Device

Each complete JAS Pediatric Knee Device requires: 1 – Frame, 1 – Thigh Cuff, 1 – Tibia Cuff



Description	Product Code
Pediatric Knee Frame	0020P-000-00

### Pediatric Knee Cuffs:

**Pediatric Knee Thigh Cuffs** (See Lower Extremity Measurement Chart on page 8 for measurement locations)

Thigh Cuff Size	Product Code		Thigh Circumference (A) in. (cm)	Thigh Length (C) in. (cm)
	Left	Right		
Pediatric	00720-P00-TL	00720-P00-TR	Up to 14 (35.6)	7 (17.9) and up
Small	00720-P0S-TL	00720-P0S-TR	14 - 17 (35.6 - 43.2)	9 3/4 (25.8) and up

**Pediatric Knee Tibia Cuffs** (See Lower Extremity Measurement Chart on page 8 for measurement locations)

Tibia Cuff Size	Product Code		Calf Circumference (B) in. (cm)	Tibia Length (D) in. (cm)
	Left	Right		
Pediatric	00720-P00-CL	00720-P00-CR	Up to 12 (30.5)	7 (17.9) and up
Xsmall	00720-PXS-CL	00720-PXS-CR	Up to 12 (30.5)	9 3/4 (25.8) and up

# LOWER EXTREMITY DEVICES

## Pediatric Ankle Device

Each complete JAS Pediatric Ankle Device requires: 1 – Frame, 1 – Calf Cuff

Description	Product Code	
Pediatric Ankle Frame	0005P-000-L1	0005P-000-R1



**Pediatric Ankle Cuffs<sup>4</sup>** (See Lower Extremity Measurement Chart on page 8 for measurement locations.)

Tibia Cuff Size	Product Code		Calf Circumference (B) in. (cm)	Tibia Length (D) in. (cm)
	Left	Right		
Small, Short	00705-PSS-CL	00705-PSS-CR	7 to 9 (17.8 – 22.9)	6 7/8 – 8 7/8 (17.5 – 22.5)
Small, Long	00705-PSL-CL	00705-PSL-CR		8 7/8 - 11 (22.5 – 27.9)
Medium, Short	00705-PMS-CL	00705-PMS-CR	9 - 11 (22.9 – 27.9)	6 7/8 – 8 7/8 (17.5 – 22.5)
Medium, Long	00705-PML-CL	00705-PML-CR		8 7/8 - 11 (22.5 – 27.9)
Large, Short	00705-PLS-CL	00705-PLS-CR	11 - 13 (27.9 – 33.0)	6 7/8 – 8 7/8 (17.5 – 22.5)
Large, Long	00705-PLL-CL	00705-PLL-CR		8 7/8 - 11 (22.5 – 27.9)
XLarge, Short	00705-PXLS-CL	00705-PXLS-CR	13 - 15 (33.0 – 38.1)	6 7/8 – 8 7/8 (17.5 – 22.5)
XLarge, Long	00705-PXLL-CL	00705-PXLL-CR		8 7/8 - 11 (22.5 – 27.9)

4.) Includes: Foot Plate Pad, Heel Cuff

# LOWER EXTREMITY DEVICES

## Lower Extremity Measurement Chart – Pediatric Devices

- A Circumference: 1 1/2" (3.8 cm) below Groin
- B Circumference: Largest point of Calf
- C Length: Groin to Knee Center (*not center of Knee Cap*)
- D Length: Knee Center (*not center of Knee Cap*) to Medial Malleolus
- E Length: Foot Length

